S300

CAT

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Table of content:

1	GENERAL INFORMATION	15
1.1	UNPACKING THE DEVICE.....	15
1.2	SAFETY INSTRUCTIONS:	15
1.2.1	<i>Explanation of the safety instructions on the device</i>	<i>17</i>
2	OPERATION AND DESCRIPTION OF THE SH10 CONTROL AND SHAKER UNIT	17
2.1	DESCRIPTION	17
2.2	PURPOSE OF THE SH10 THERMO SHAKER.....	17
2.3	DIMENSIONS AND INSTALLATION OF THE SH10 SHAKER	18
2.3.1	<i>Ventilation requirements during operation</i>	<i>18</i>
2.4	OPERATION OF THE S300 SHAKER	19
2.5	MALFUNCTIONS OF THE DEVICE	19
3	MAINTENANCE AND SERVICE	20
3.1	SERVICE UND REPAIR.....	20
3.2	CLEANING AND DECONTAMINATION OF THE DEVICE	20
4	DE-INSTALLATION AND DISPOSAL OF THE SH10 THERMO SHAKER	20
4.1	DE-INSTALLATION OF THE DEVICE.....	21
4.2	DISPOSAL OF THE DEVICE	21
5	PACKING AND STORAGE OF THE DEVICE	21
5.1	PACKING THE DEVICE	21
5.2	STORING THE DEVICE	21
6	WARRANTY AND EXCLUSION OF LIABILITY	22
7	TECHNICAL DATA	23

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1 General Information

1.1 Unpacking the device

Check the outside box around the instrument. If there is any damage such as pushed-in sides, chemical damage to the box, water marks or other physical effects that may have damaged the contents, then inform the forwarding agent and ask about the necessary measures to be taken. If there is no visible damage on the outside of the box then open it carefully.

Please check the contents for damage or deviation from the packing list:

1x S300 (Artikel-Nr.: 60388-00/ 60388-01)
1x Operating Manual

Note:

If all parts are complete then you can start to operate the device after reading the instructions.

Note:

Please read the operating instructions thoroughly and make sure that you have understood everything before starting the device. Should you have any questions about installation, operation or maintenance after reading the instructions then please contact the supplier or manufacturer of the device at the following address:

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1.2 Safety instructions:

1. Please read the user manual completely before starting the device. Points, that require particular attention for your application should be marked so that they can be easily seen.
2. These instructions should always be at hand particularly for those who want to use the device. Hence the instructions should be kept close to the device.
3. Only use original parts and accessories from the manufacturer. These parts are matched to the tight tolerances of the devices and ensure reliability.

The characteristics of the device and especially the use of selected materials, ensure versatile use in the lab. Details about the resistance of the device against external mediums can be found in the section: "Technical data".

Explanation of the safety instructions in the operating manual:

	This symbol highlights the following instructions, which the user has to follow strictly! Failure to follow these instructions can impair the safe function and safety of the user.
	This symbol highlights bans, which the user has to follow strictly! Failure to follow these bans can seriously impair the functionality and safety of the user.
	This symbol highlights instructions that should be strictly followed by the user to ensure safe operation of the device.
	This warning symbol indicates that the heating surface is not to be touched during operation. The hotplate of the SH10 shaker unit can reach 95°C during operation. This means that there is a risk of burning for the user.

Note:

These instructions don't cover all safety risks that may result from using the device or from chemicals, reagents, devices and equipment that are used in a specific test or process. It is the users responsibility to get advice and stipulate the correct health and safety regulations or practices and also the use or limitations of this device before using it.

Note:

Failure to follow the safety and operation instructions in this manual and all incorrect handling of the device, can impair the foreseen protection of the device as well as void the warranty!

Caution:

Only instructed users may start the device.

Caution:

Don't use this device in dangerous surroundings and especially not in explosive ones.

Caution:

The main switch of the control box does not fully disconnect the device from the power source. Only one phase is separated. Please pull out the mains plug to disconnect the device completely from the mains/ power supply.

Note:

Make sure that all safety instructions of the device or on the device can be clearly seen during operation.

Caution:

Don't open the device. Repairs are only to be carried out by trained service technicians.

Caution:

The shaker may not be connected to or disconnected from the control unit as long as the main switch of the control unit is on. This can damage the device! We recommend to connect the shaker with the control unit before connecting the power supply.

Caution:

Don't use flammable media together with or in the direct surroundings of the device.

1.2.1 Explanation of the safety instructions on the device

Warning against general dangers: This symbol indicates that it is imperative to read and understand the operating instructions of the device before starting it. Incorrect handling can endanger the safe operation of the device and health of the user.

2 Operation and description of the SH10 control and shaker unit

2.1 Description

The S300 is a shaking device with an orbital shaking movement. The device will be powered on and off by switching the power-button located at the frontside of the housing. The shaking frequency is factory-provided fixed to 300 rpm. The shaking frequency cannot be changed by the user.

2.2 Purpose of the SH10 thermo shaker

The SH10 thermo shaker is designed for the following tasks. All incorrect handling can lead to malfunction of the device, or even endanger the user.

The S300 shaker is for mixing/ shaking of media. These media may only be placed on the shaker using suitable reagent containers. The placing of media directly on the platform is strictly forbidden. The reagent or medium containers may not be flammable and have to be fixed on the shaker platform. Ensure that the reagent or medium containers are suitable for holding your medium.

Caution:

Don't use flammable media together with or in the direct surroundings of the device.

Caution:

It is the user task to find out whether the device is suitable for the application. If in doubt clarify this with your dealer or directly with the manufacturer.

Caution:

Only use the device so that there is no danger to the user or anyone else. Avoid splashes use safety glasses, protective clothing and safety gloves.

2.3 Dimensions and installation of the SH10 shaker

The installation dimensions such as length x width and height, are shown in the following figure:

Caution:

It is essential to install the shaker module on a suitable base plate or holder, that is secured against slipping. Don't start the device as long as it is not adequately positioned.

Caution:

During installation make sure that the device is at a safe distance from its direct surroundings. The shaker unit should be at a safe distance of least 2 cm from its surroundings. The upward safety distance should be 30 cm due to the hotplate. The mentioned safety distances only refer to the functionality of the device and not to any dangers in relation to fire protection. Don't use any flammable media right beside the device.

Caution:

During installation of the device make sure that the power and control cables are correctly laid. No lines may run directly beside the hotplate. Keep a safe distance of 2 cm at the side in this case. Don't lay any cables on the heating surface.

Note:

When installing the device and power supply make sure that the associated circuit breakers are freely accessible!

2.3.1 Ventilation requirements during operation

Ensure adequate and continuous ventilation of the SH10 thermo shaker during operation. The ambient temperature of the shaker may not exceed 40°C. You as the user have to take care that the device does not overheat by providing adequate ventilation.

The user has to provide the ventilation. Whether the ventilation has an effect on the accuracy of the temperature and temperature stability of the shaker unit, and whether these have a negative effect on the method/ application has to be judged by the user.

Failure to follow the instructions means that the safe operation of the device cannot be guaranteed and that there may be a danger to people and material!

2.4 Operation of the S300 shaker

Make sure that the device is only used in the ambient conditions as described in the “Technical Data” section. Please read the operating instructions before starting the device. Please follow all instructions mentioned on the operating instructions.

As soon as the device is connected to an adequate power source, it can be powered on. By switching on the power-button at the frontside of the device, it starts shaking the platform. By switching off the power-button, the shaking platform will stop moving at once.

Caution:

The main switch of the control unit does not fully disconnect the device from the power source when switched off. The main switch only disconnects one phase. Please pull out the plug to disconnect the device completely from the mains.

Achtung:

Please do not use any other power-source than stated on the product label on the device. These devices can be purchased as 110V/60Hz or as 230V/50Hz version.

Caution:

It is the user task to find out whether the device is suitable for the application. If in doubt clarify this with your dealer or directly with the manufacturer.

Caution:

The placing of media directly on the shaking platform is strictly forbidden. The reagent or medium containers may not be flammable and have to be fixed on the shaker platform. Ensure that the reagent or medium containers are suitable for holding your medium.

2.5 Malfunctions of the device

Function	Potential error	Possible reasons	Possible effects	Possible checks of the individual functions	Possible countermeasures
Shaking of the platform	Device does not shake	1) Device is not connected to mains 2) Shaking platform is blocked 3) Load on shaking platform is too high	No known effects on the device. User to determine the effects this may have on the application	Visible status query (by eye or by external rpm-measuring device)	Switch the device off correctly for a brief period and then on again via the mains switch of the control unit. Before doing so make sure that no filled reagent container is on the shaker platform. Protect yourself from splashing of / or direct contact with the medium you are using.
	Device shakes uncontrollably or with the wrong speed	1) Shaking platform is blocked 2) Faulty electronic	No known effects on the device. User to determine the effects this may have on the application	Visible status query (by eye or by external rpm-measuring device)	If the device still does not work correctly then contact the manufacturer or supplier of the device.

Caution:

Don't open the device. Repairs are only to be carried out by trained service technicians. Please contact the manufacturer or supplier for repair or maintenance.

3 Maintenance and service

3.1 Service und Repair

Should you have any questions about installation, operation or maintenance then please contact the supplier or manufacturer of the device at the following address:

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E-Mail: info@cat-ing.de

Web: www.cat-ing.de

Don't carry out any electronic or mechanical repairs on the device. There are no parts inside the device that can be maintained or repaired by the user. Repairs by the user may damage the electronics. The warranty may also be void as a result.

Caution:

Don't open the device! Repairs are only to be carried out by trained service technicians.

Caution:

In case of repair: If you are obliged to send the device back to the manufacturer, please make sure that the device has been adequately cleaned / decontaminated beforehand. The user is responsible for direct or indirect damage both corporal and material, that may result from the media used.

3.2 Cleaning and decontamination off the device

Wipe external surfaces and controls of the device with a mild soap solution (water and household detergent) and a lint-free cloth. Just use a small amount of cleaning agent and not wet cloths. Under no circumstances should the surfaces be cleaned with acids, chlorine containing solvents or saline solutions.

Please don't scour or scratch, otherwise damage to the lettering and surfaces cannot be avoided. Should you wish to use cleaners other than those mentioned make sure beforehand that the foreseen method will not cause any damage.

4 De-installation and disposal of the SH10 thermo shaker

Correct handling is required due to the device characteristics and the possible area of application. Please observe the following instructions when using the device:

4.1 De-installation of the device

When de-installing / disassembling the device and its components it is essential to take the following steps:

1. First switch off the S300 device with the mains switch.
2. Disconnect the device from the mains
5. Now you can remove the shaker and its components from the work area.

4.2 Disposal of the device

Damaged or disposed electric or electronic devices have to be handed in at recycling centres. The packing material too should be disposed of properly (material separation).

Alternatively we can carry out proper disposal of the device for you. Please contact us for this!

5 Packing and storage of the device

5.1 Packing the device

Caution:

Before packing make sure that the device has been adequately cleaned / decontaminated. User who unpack the device afterwards and want to start it up should be aware of potential dangers. The user is responsible for direct or indirect damage both corporal and material, that may result from the media used as long as they don't provide adequate signs about the potential dangers.

When possible use the original material for packing the device. Protect each module against shocks from outside with bubble wrap. If you pack individual modules in a box ensure adequate spacing between them and fill the spaces with insulation material. If you don't use the original packing then mark the pack on the outside as follows:

- Glass sign (Caution, fragile)
- Umbrella (keep dry)
- Pack content (details about the content)

5.2 Storing the device

Observe the following notes for correct storage of the thermo shaker.

1. Separate the individual modules from each other for storing.
2. Pack each module in a separate airtight plastic bag.

3. Ensure the following ambient conditions:

Max. ambient temperature: 5-40°C
Max. air humidity: 80%

Caution:

Before packing make sure that the device has been adequately cleaned / decontaminated.

User who unpack the device afterwards and want to start it up should be aware of potential dangers. The user is responsible for direct or indirect damage both corporal and material, that may result from the media used as long as they don't provide adequate signs about the potential dangers.

6 Warranty and exclusion of liability

The manufacturer commits to repair or replace this product, at his discretion, in case of defects arising from material or workmanship, that appear within 24 months of delivery of said product. In case of replacement the warranty continues until the 24 months have elapsed from the date of purchase.

The dealer or manufacturer is to be involved for all customer service matters.

This warranty is not valid in case of defects or malfunctions caused by accidents, negligence, incorrect use, wrong servicing and other grounds that cannot be traced back to faulty materials or workmanship.

All warranty and guarantee rights are limited to correct and proper use of the device while following the general safety regulations and these operating instructions. The manufacturer is not obliged to provide other performances apart from those contained in this warranty declaration. In particular corporal or material damages resulting from defective parts or other functional errors are excluded from liability.

The manufacturer's obligation is limited in scope to the repair or replacement of faulty parts. Under no circumstances is the manufacturer obliged to pay for all types of damages resulting from the use and ownership of this product.

Caution:

Please don't try to repair the device in case of a breakdown. There are no parts inside the device that can be maintained or repaired by the user. Repair by the user can void the warranty.

Caution:

Don't open the device. Repairs are only to be carried out by trained service technicians.

7 Technical Data

SH10 Controller	65581-00
Power Supply	24V DC +/-10% 90W
Motor Control Speed	<ul style="list-style-type: none"> – continuously adjustable – microprocessor controlled – Soft-Start of motor
Ambient temperature range	5 - 40 °C
Max. humidity	80 % relative humidity
Dimensions (LxWxH)	170mm x 140mm x 90mm
Weight	1,6 Kg
Pollution Degree	II
Remark	For Indoor Use Only! 

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